ZENET

MASSAGE CUSHION ZET-759

Massage Cushion ZET-759

INSTRUCTIONS FOR USE

Power 45W

Voltage 220V

10 min timer

Dear customer,

Thank you for choosing the TM ZENET device

Please read the operating instructions and keep them for future reference.

Attention

Remove the packaging and make sure that the device is not damaged. If in doubt, do not use the device and contact a service center. Keep packaging materials such as bags and staples out of the reach of children. Keep children away from using the massager as a toy.

The massage device can only be repaired by a qualified specialist. Do not try to disassemble the massager for details. Any contact with water is prohibited (except for cleaning with a damp cloth). Never immerse the device in water, do not use it in the bathroom, in the shower or in the pool.

Consult a doctor before using the massager:

- during pregnancy
- if you have a pacemaker
- if the following diseases are present (e.g. varicose veins, thrombosis, jaundice, diabetes, sciatica, intervertebral hernia, etc.)
- during the recovery phase after the operation.

Persons (including children) with limited physical or mental abilities, as well as persons who do not have the experience or knowledge required to operate this massager, may only use the massager if they are under the supervision of persons who are responsible for their safety.

Switch off the massager immediately if you experience pain or discomfort.

If you experience pain in the muscles or joints after a massage, refuse to use the massager and consult a doctor.

Attention

The massager is not a medical device. Massager is not a substitute for medical

or physiotherapy treatment. If you have any doubts about the possibility of using a massager, consult a doctor.

FUNCTIONS:

1. Plug in the device. Press the "On / Off" button to start the massager. By default, the device is switched on in the soothing massage mode.
2. Press the warm-up button to activate the warm-up function. The heating is switched off by pressing it again.
3. When the mode selection button is pressed, the device switches to the pulsed massage mode, and when it is pressed again, the massager returns to the calming massage mode.
4. With the Auto button you can start one of the automatic massage programs. The first press starts the Deep Relaxing Massage program, the second press starts the Relaxation Massage program.
5. At the end of the massage, switch off the massager with the "On / Off" button and pull the mains plug out of the socket.

Cleaning

Remove the plug from the socket for cleaning. After cooling, wipe the massager with a slightly damp cloth. Make sure that water never gets into the device. Do not use the device until it is completely dry.

GUARANTEE LIABILITIES:

1. Guarantee of the device according to the operating rules - 1 year from the date of sale.
2. During the warranty period, the device must be replaced if any defects caused by the manufacturer are found. The sales date and the stamp of the trade organization are required.
3. The guarantee does not apply to devices with:
 - mechanical damage;
 - defects resulting from a violation of the terms of use or incorrect actions by the owner of the device;
 - deficiencies due to natural disasters;
 - defects caused by the ingress of foreign bodies and liquids as a result of self-repair or self-introduction of structural changes;
 - used for commercial or industrial purposes.

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